

BRUNO VASSARI

LAB DIVISION

HA50X PRO-HYALURONIC TREATMENT

TREATMENT WITH HYALURONIC ACID



Hyaluronic acid is one of the main components of our skin.

As the years go by the amount that our organism synthesises decreases and the skin is affected; dehydration, wrinkles, lack of tone and density are some of these effects.

This treatment aims to boost the levels of hyaluronic acid to maintain the skin in its best condition.

INDICATIONS OF THE TREATMENT:

Indicated for all types of skin, especially those that need to re-establish their degree of moisturisation.

FREQUENCY OF TREATMENT:

Undertake a weekly session for 6 weeks.

STEP BY STEP

(Duration of treatment: 60 minutes)

- 1 We start the treatment by removing the make-up from the eyes and lips gently with cotton wool soaked in **Bio Douce**.
- 2 We then cleanse and tone up with the best products adapted to the type of skin.
- 3 With the skin now cleansed, we apply a generous layer of **Hyaluronic Scrub and Mask** avoiding the contour of the eyes. Once applied, we make gentle circular movements dampening the fingertips repeatedly to facilitate movement. After the exfoliation, we leave the product exposed for a few minutes and remove it using sponges and abundant warm water.
- 4 Then we give a facial massage with **Vitalizing Massage Cream** applying the cream with a fan brush and leaving it exposed for a few minutes. Then we give a facial massage for 15-20 minutes.
- 5 After completing the massage, we remove the cream with a towel dampened in hot water.
- 6 Now with the skin clean, we apply **Hyaluronic Serum** and massage it until it is completely absorbed.
- 7 We then apply the **Pro-Hyaluronic Mask**. To prepare it we pour 70-80 ml of cold water into a bowl and sprinkle the contents of the sachet over it. Leave to rest until the powder penetrates the water and then we mix it until obtaining a mousse texture that we apply onto the face immediately, avoiding the contour of the eyes. Leave exposed for 15 minutes.
- 8 After this exposure time, we remove the excess of the product with a spatula and we complete the cleansing with some mitts dampened in warm water.
- 9 We finish the treatment with the **Hyaluronic Cream Pollution Neutralizer** finalising cream.

Follow-up at home:

Exfoliate

Hyaluronic Scrub and Mask

Double-action mask that enables you to exfoliate and treat the skin. Due to the combination of exfoliant particles of perlite and white clay the skin remains clean, mattified and free of dead cells.

Apply a generous layer and massage dampening the fingertips. Leave exposed for a few minutes and remove with abundant warm water.

Day treatment

Hyaluronic Cream Pollution Neutralizer

Formulated anti-pollution cream, both for the daily protection of the skin against harmful agents and for strengthening its natural defence mechanisms. It includes two types of hyaluronic acid, of low and high molecular weight, which ensures the maximum moisturising of the skin, its turgidity and flexibility.

Apply in the morning after regular cleansing.

Intensive treatment

Multi molecular Hyaluronic Filler

Extra-concentrated hyaluronic acid serum that provides intense moisturising throughout the day. The skin feels more rejuvenated and firm, elastic and flexible.

Apply morning and/or night before the regular cream treatment.

Weekly treatment

Hyaluronic Sleeping Mask

Night mask formulated with hyaluronic acid and a polysaccharide that treats the skin during sleep to wake up with a rested, nourished and more uniform face. It improves the moisturisation of the skin and protects it from external aggressions.

Apply a fine layer before going to sleep. Apply daily as an intensive treatment.